



M18	1	DIA42X2 TORIC SEAL	JOINT FRANCAIS	133107
	1			
M16	1	MALE SMS51 SANITARY COUPLING	316L St STEEL	-
M15	1			-
M14	1	FEMALE SMS51 SANITARY COUPLING	316L St STEEL	-
-	1	-	-	-
-	-	-	-	-
M10	8	FLAT SLICE	304L St STEEL	
M9	8	HEXAGONAL M12 BOLT	304L St STEEL	
M8	8	GEARBOX ASSEMBLY BOLT	304L St STEEL	SCREW HC 12 x 4
M7	1	SHAFT ASSEMBLY SCREW	316L St STEEL	HEAD H 8 x 30
-	-	-	-	-
-	-	-	-	-
M4	1	PADDLE SEAL	BAUM	75 x 95 x 12
M3	1	SHAFT SEAL	BAUMSL	38 x 55 x 7
M2	1	SHAFT PROTECTION	BONFIGLIOLI	
M1	1	MOTOREDUCER 1.1 KW 30Tr/Mn	BONFIGLIOLI	
8	1	FLANGE SEAL	WHITE RUBBER	070619-100-17
7	1	FLAT SLICE	INOX 316L	DP-400-100-16
6	1	ASSEMBLY FLANGE	INOX 316L	DP-400-100-15
5	1	COVER	INOX 316L	DP-400-100-14
4	1	PADDLE ASSEMBLY SCREW	INOX 316L	DP-400-100-13
3	1	SHAFT	INOX 316L	DP-400-100-12
2	1	PADDLE	INOX 316L	DP-400-100-11
1	1	DP400 BODY	INOX 316L	DP 400-100-10
N°	Qt?	D?signation	Observation	Plan

TECHNI lab		CE2948	F91029	EVERY Cedex
CE PLAN EST LA PROPRIETE DE TECHNI lab ET A CE TITRE IL NE DOIT ETRE COPIE, TRANSMIS A UNE TIERCE PERSONNE OU MODIFIE SANS AUTORISATION ECRITE DE TECHNI lab		ATTENTION TOUTES LES DIMENSIONS SONT INDICATIVES ET PEUVENT CHANGER SANS PREAVIS. NE PAS UTILISER POUR LA FABRICATION		
Client :				
Projet : -				
Titre ENSEMBLE_DP_400_BAS				INDICE: F
Plan N° : DP-400-100-00	Dessin? par :	Format :		
Le: 14/09/07	Ech.: /	M.BRAS	A3	

